

OVATION™ HIP STEM

PRODUCT SUMMARY




Ortho Development introduces the Ovation™ Hip Stem, an innovation in design and simplicity that reduces bone loss and length of surgery. To provide options for various patient needs, the Ovation is available in a range of sizes with both standard and extended neck offset options.

- Tapered stem design provides excellent rotational stability and promotes 3-point fixation.
- Minimal instrumentation simplifies surgical procedure.
- Streamlined profile is bone conserving.
- Optimized neck geometry with highly polished surface improves overall range of motion and decreases wear.
- Available in 12 sizes with standard and extended neck offsets to accommodate patient anatomy.
- Proximal body surface is coated with titanium plasma spray for solid, long-term fixation.
- Distal body flute to more closely match bone properties.
- Accommodates all surgical approaches, including direct anterior approach.

PRODUCT SPECIFICATIONS

SIZE	STEM LENGTH	PROXIMAL WIDTH (M/L)	DISTAL WIDTH	NECK LENGTH	HORIZONTAL OFFSET (STANDARD)	HORIZONTAL OFFSET (EXTENDED)
1	98	26	5	29	34	39
2	106	27	6	30	34	39
3	108	28	6	30	35	40
4	110	29	7	32	37	43
5	114	30	8	32	38	44
6	118	31	9	33	38	44
7	122	32	10	34	42	49
8	126	34	11	34	43	50
9	130	36	13	34	44	51
10	134	38	16	36	47	55
11	138	40	18	36	48	56
12	142	42	20	36	49	57

Ovation Stem Dimensions and Offsets (all dimensions are in mm)



ORDER INFORMATION

O V A T I O N H I P S T E M

PRODUCT NO.	DESCRIPTION	PRODUCT NO.	DESCRIPTION
101-0001	Ovation Stem Standard Size 1	101-1001	Ovation Stem Extended Size 1
101-0002	Ovation Stem Standard Size 2	101-1002	Ovation Stem Extended Size 2
101-0003	Ovation Stem Standard Size 3	101-1003	Ovation Stem Extended Size 3
101-0004	Ovation Stem Standard Size 4	101-1004	Ovation Stem Extended Size 4
101-0005	Ovation Stem Standard Size 5	101-1005	Ovation Stem Extended Size 5
101-0006	Ovation Stem Standard Size 6	101-1006	Ovation Stem Extended Size 6
101-0007	Ovation Stem Standard Size 7	101-1007	Ovation Stem Extended Size 7
101-0008	Ovation Stem Standard Size 8	101-1008	Ovation Stem Extended Size 8
101-0009	Ovation Stem Standard Size 9	101-1009	Ovation Stem Extended Size 9
101-0010	Ovation Stem Standard Size 10	101-1010	Ovation Stem Extended Size 10
101-0011	Ovation Stem Standard Size 11	101-1011	Ovation Stem Extended Size 11
101-0012	Ovation Stem Standard Size 12	101-1012	Ovation Stem Extended Size 12



12187 South Business Park Drive
Draper, Utah 84020

Telephone: 801.553.9991
Fax: 801.553.9993

For more information visit
www.orthodevelopment.com