

PIVOT BIPOLAR™



O R T H O D E V E L O P M E N T


INTRODUCTION TO PIVOT BIPOLAR™

The Pivot Bipolar™ was designed with a simple one-step assembly interface between the femoral head and the compression molded liner which requires no instruments and reduces OR time.


FEATURES & BENEFITS

- Compression Molded UHMWPE
- Available in 38-60mm (1mm increments)
- Simple one-step assembly
- 75° Range of Motion


ASSEMBLY AND PULL-OUT STRENGTH



PRODUCT SPECIFICATIONS



DESIGN RATIONALE



Once the appropriate femoral head is placed onto the stem the Pivot Bipolar is then placed onto the femoral head. The Pivot Bipolar is then advanced until it is past the retaining ring and snapped into place.

RANGE OF MOTION

75° Range of Motion
(12/14 taper and 28mm femoral head)

ORDER INFORMATION

SIZE	OUTER DIA.	HEAD DIA.	SIZE	OUTER DIA.	HEAD DIA.
22x38mm	38mm	22mm	22x39mm	39mm	22mm
22x40mm	40mm	22mm	28x41mm	41mm	28mm
28x42mm	42mm	28mm	28x43mm	43mm	28mm
28x44mm	44mm	28mm	22x45mm	45mm	22mm
28x46mm	46mm	28mm	28x47mm	47mm	28mm
28x48mm	48mm	28mm	28x49mm	49mm	28mm
28x50mm	50mm	28mm	28x51mm	51mm	28mm
28x52mm	52mm	28mm	28x53mm	53mm	28mm
28x54mm	54mm	28mm	28x55mm	55mm	28mm
28x56mm	56mm	28mm	28x57mm	57mm	28mm
28x58mm	58mm	28mm	28x59mm	59mm	28mm
28x60mm	60mm	28mm			

HEAD OD	OFFSET	HEAD OD	OFFSET	HEAD OD	OFFSET
22mm	-3mm	22mm	+0mm	22mm	+3mm
22mm	+6mm	28mm	-3mm	28mm	+0mm
28mm	+3mm	28mm	+6mm	28mm	+9mm
32mm	-3mm	32mm	+0mm	32mm	+3mm
32mm	+6mm	32mm	+9mm		

PRODUCT NO.	DESCRIPTION	PRODUCT NO.	DESCRIPTION
133-2238	Pivot Bipolar Cup 22x38mm	133-2850	Pivot Bipolar Cup 28x50mm
133-2239	Pivot Bipolar Cup 22x39mm	133-2851	Pivot Bipolar Cup 28x51mm
133-2240	Pivot Bipolar Cup 22x40mm	133-2852	Pivot Bipolar Cup 28x52mm
133-2241	Pivot Bipolar Cup 22x41mm	133-2853	Pivot Bipolar Cup 28x53mm
133-2242	Pivot Bipolar Cup 22x42mm	133-2854	Pivot Bipolar Cup 28x54mm
133-2243	Pivot Bipolar Cup 22x43mm	133-2855	Pivot Bipolar Cup 28x55mm
133-2244	Pivot Bipolar Cup 22x44mm	133-2856	Pivot Bipolar Cup 28x56mm
133-2245	Pivot Bipolar Cup 22x45mm	133-2857	Pivot Bipolar Cup 28x57mm
133-2846	Pivot Bipolar Cup 28x46mm	133-2858	Pivot Bipolar Cup 28x58mm
133-2847	Pivot Bipolar Cup 28x47mm	133-2859	Pivot Bipolar Cup 28x59mm
133-2848	Pivot Bipolar Cup 28x48mm	133-2860	Pivot Bipolar Cup 28x60mm
133-2849	Pivot Bipolar Cup 28x49mm		

PRODUCT NO.	DESCRIPTION	PRODUCT NO.	DESCRIPTION
138-2230	CoCr Femoral Head 22mm -3	138-2806	CoCr Femoral Head 28mm +6
138-2200	CoCr Femoral Head 22mm +0	138-2809	CoCr Femoral Head 28mm +9
138-2203	CoCr Femoral Head 22mm +3	138-3230	CoCr Femoral Head 32mm -3
138-2206	CoCr Femoral Head 22mm +6	138-3200	CoCr Femoral Head 32mm +0
138-2830	CoCr Femoral Head 28mm -3	138-3203	CoCr Femoral Head 32mm +3
138-2800	CoCr Femoral Head 28mm +0	138-3206	CoCr Femoral Head 32mm +6
138-2803	CoCr Femoral Head 28mm +3	138-3209	CoCr Femoral Head 32mm +9



12187 South Business Park Drive
Draper, Utah 84020

Telephone: 801.553.9991
Fax: 801.553.9993

For more information visit
www.orthodevelopment.com