

SOLITUDE™
UNIPOLAR




O R T H O D E V E L O P M E N T

INTRODUCTION TO SOLITUDE™ UNIPOLAR

The Solitude™ Unipolar features a simple intraoperative assembly. It is available in sizes 38-60mm (1mm increments). Also the trunion adapters are available in -3, +0, +3, +6, and +9mm options.

PRODUCT SPECIFICATIONS



See back panel for sizing.



SIZE	OUTER DIAMETER	SIZE	OUTER DIAMETER
22x38mm	38mm	22x39mm	39mm
22x40mm	40mm	28x41mm	41mm
28x42mm	42mm	28x43mm	43mm
28x44mm	44mm	22x45mm	45mm
28x46mm	46mm	28x47mm	47mm
28x48mm	48mm	28x49mm	49mm
28x50mm	50mm	28x51mm	51mm
28x52mm	52mm	28x53mm	53mm
28x54mm	54mm	28x55mm	55mm
28x56mm	56mm	28x57mm	57mm
28x58mm	58mm	28x59mm	59mm
28x60mm	60mm		

12/14 TRUNION ADAPTER OFFSET
- 3mm
+0mm
+3mm
+6mm
+9mm

ORDER INFORMATION

PRODUCT NO.	DESCRIPTION	PRODUCT NO.	DESCRIPTION
132-0038	Solitude Unipolar Head 38mm	132-0050	Solitude Unipolar Head 50mm
132-0039	Solitude Unipolar Head 39mm	132-0051	Solitude Unipolar Head 51mm
132-0040	Solitude Unipolar Head 40mm	132-0052	Solitude Unipolar Head 52mm
132-0041	Solitude Unipolar Head 41mm	132-0053	Solitude Unipolar Head 53mm
132-0042	Solitude Unipolar Head 42mm	132-0054	Solitude Unipolar Head 54mm
132-0043	Solitude Unipolar Head 43mm	132-0055	Solitude Unipolar Head 55mm
132-0044	Solitude Unipolar Head 44mm	132-0056	Solitude Unipolar Head 56mm
132-0045	Solitude Unipolar Head 45mm	132-0057	Solitude Unipolar Head 57mm
132-0046	Solitude Unipolar Head 46mm	132-0058	Solitude Unipolar Head 58mm
132-0047	Solitude Unipolar Head 47mm	132-0059	Solitude Unipolar Head 59mm
132-0048	Solitude Unipolar Head 48mm	132-0060	Solitude Unipolar Head 60mm
132-0049	Solitude Unipolar Head 49mm		

PRODUCT NO.	DESCRIPTION	PRODUCT NO.	DESCRIPTION
132-123012/14	Solitude Unipolar Adapter - 3	132-120612/14	Solitude Unipolar Adapter +6
132-120012/14	Solitude Unipolar Adapter +0	132-120912/14	Solitude Unipolar Adapter +9
132-120312/14	Solitude Unipolar Adapter +3		



12187 South Business Park Drive
Draper, Utah 84020

Telephone: 801.553.9991
Fax: 801.553.9993

For more information visit
www.orthodevelopment.com