

## MODULAR TIBIAL COMPONENTS

- Available in sizes 1-6
- Accommodates both primary PS and revision SCPS\* tibial inserts
- Available in 0° neutral slope and 5° posterior slope
- One tibial keel length simplifies reaming process
- Patented Balanced Knee System tibial insert locking mechanism\*\*
- Solid tibial base eliminates pathways for debris migration
- Augments available in hemi wedges and hemi blocks
- Lateral tapered edge on augments matches tibial anatomy
- Uses same modular stems as modular femoral components
- Stems available in both cemented and fluted options
- SnapLoc™ feature on augments for simple initial fixation
- Spirallock® threads on stems and augments for secure fixation



## SEMI-CONSTRAINED POSTERIOR STABILIZED (SCPS) TIBIAL INSERT

- Available in sizes 1-6 in thickness variations from 8-30mm
- Patented Balanced Knee System tibial insert locking mechanism\*\*
- SCPS post designed with central pin for strength
- Patellar notch designed for increased patellar relief



\*SCPS tibial insert may only be used with modular tibial tray  
\*\*US Patent No. 6,379,388

## ORDER INFORMATION (CONTINUED)

### TIBIAL AUGMENTS (CONTINUED)

PRODUCT NO.	DESCRIPTION	PRODUCT NO.	DESCRIPTION
567-6201	Tibial Augment 10mm Block Lt Lat/Rt Med Size 1	567-7104	Tibial Augment 15mm Block Rt Lat/Lt Med Size 4
567-6202	Tibial Augment 10mm Block Lt Lat/Rt Med Size 2	567-7105	Tibial Augment 15mm Block Rt Lat/Lt Med Size 5
567-6203	Tibial Augment 10mm Block Lt Lat/Rt Med Size 3	567-7106	Tibial Augment 15mm Block Rt Lat/Lt Med Size 6
567-6204	Tibial Augment 10mm Block Lt Lat/Rt Med Size 4	567-7201	Tibial Augment 15mm Block Lt Lat/Rt Med Size 1
567-6205	Tibial Augment 10mm Block Lt Lat/Rt Med Size 5	567-7202	Tibial Augment 15mm Block Lt Lat/Rt Med Size 2
567-6206	Tibial Augment 10mm Block Lt Lat/Rt Med Size 6	567-7203	Tibial Augment 15mm Block Lt Lat/Rt Med Size 3
567-7101	Tibial Augment 15mm Block Rt Lat/Lt Med Size 1	567-7204	Tibial Augment 15mm Block Lt Lat/Rt Med Size 4
567-7102	Tibial Augment 15mm Block Rt Lat/Lt Med Size 2	567-7205	Tibial Augment 15mm Block Lt Lat/Rt Med Size 5
567-7103	Tibial Augment 15mm Block Rt Lat/Lt Med Size 3	567-7206	Tibial Augment 15mm Block Lt Lat/Rt Med Size 6

### FEMORAL AUGMENTS

PRODUCT NO.	DESCRIPTION	PRODUCT NO.	DESCRIPTION
568-1104	Femoral Augment Distal Size 1 4mm	568-1704	Femoral Augment Distal Size 7 4mm
568-1108	Femoral Augment Distal Size 1 8mm	568-1708	Femoral Augment Distal Size 7 8mm
568-1204	Femoral Augment Distal Size 2 4mm	568-1712	Femoral Augment Distal Size 7 12mm
568-1208	Femoral Augment Distal Size 2 8mm	568-2104	Femoral Augment Posterior Size 1 4mm
568-1212	Femoral Augment Distal Size 2 12mm	568-2108	Femoral Augment Posterior Size 1 8mm
568-1304	Femoral Augment Distal Size 3 4mm	568-2204	Femoral Augment Posterior Size 2 4mm
568-1308	Femoral Augment Distal Size 3 8mm	568-2208	Femoral Augment Posterior Size 2 8mm
568-1312	Femoral Augment Distal Size 3 12mm	568-2304	Femoral Augment Posterior Size 3 4mm
568-1404	Femoral Augment Distal Size 4 4mm	568-2308	Femoral Augment Posterior Size 3 8mm
568-1408	Femoral Augment Distal Size 4 8mm	568-2404	Femoral Augment Posterior Size 4 4mm
568-1412	Femoral Augment Distal Size 4 12mm	568-2408	Femoral Augment Posterior Size 4 8mm
568-1504	Femoral Augment Distal Size 5 4mm	568-2504	Femoral Augment Posterior Size 5 4mm
568-1508	Femoral Augment Distal Size 5 8mm	568-2508	Femoral Augment Posterior Size 5 8mm
568-1512	Femoral Augment Distal Size 5 12mm	568-2604	Femoral Augment Posterior Size 6 4mm
568-1604	Femoral Augment Distal Size 6 4mm	568-2608	Femoral Augment Posterior Size 6 8mm
568-1608	Femoral Augment Distal Size 6 8mm	568-2704	Femoral Augment Posterior Size 7 4mm
568-1612	Femoral Augment Distal Size 6 12mm	568-2708	Femoral Augment Posterior Size 7 8mm

### JUNCTION BOX

PRODUCT NO.	DESCRIPTION
561-5000	Femoral Junction Box 5°



12187 South Business Park Drive  
Draper, Utah 84020

Telephone: 801.553.9991  
Fax: 801.553.9993

For more information visit  
[www.orthodevelopment.com](http://www.orthodevelopment.com)

Ortho Development and all other registered trademarks and trademarks are the property of Ortho Development Corporation. Patent(s) pending.  
Part No. 1030 REV 09/08

## BALANCED KNEE® SYSTEM REVISION

The Balanced Knee System Revision utilizes the same patented locking mechanism and balancing techniques as the primary Balanced Knee System. The Balanced Knee System Revision is designed to provide surgeons with several options in both primary and revision situations. This includes an easy intraoperative transition from the primary Balanced Knee System components to Balanced Knee System Revision if necessary.

## MODULAR FEMORAL COMPONENTS

- Available in sizes 1-7
- Articulates on both primary Posterior Stabilized (PS) and revision Semi-constrained Posterior Stabilized (SCPS®) tibial inserts
- Uses same modular stems as modular tibial components
- Stems available in both cemented and fluted options
- Distal and posterior augments available
- SnapLoc™ feature on augments for simple initial fixation
- Spirallock® threads on stems and augments for secure fixation
- Junction box allows for A/P translation of stem



## ORDER INFORMATION

### MODULAR FEMORAL COMPONENTS

PRODUCT NO.	DESCRIPTION	PRODUCT NO.	DESCRIPTION
561-1101	Modular Femoral Component Left Size 1	561-1102	Modular Femoral Component Right Size 1
561-1201	Modular Femoral Component Left Size 2	561-1202	Modular Femoral Component Right Size 2
561-1301	Modular Femoral Component Left Size 3	561-1302	Modular Femoral Component Right Size 3
561-1401	Modular Femoral Component Left Size 4	561-1402	Modular Femoral Component Right Size 4
561-1501	Modular Femoral Component Left Size 5	561-1502	Modular Femoral Component Right Size 5
561-1601	Modular Femoral Component Left Size 6	561-1602	Modular Femoral Component Right Size 6
561-1701	Modular Femoral Component Left Size 7	561-1702	Modular Femoral Component Right Size 7

### MODULAR TIBIAL TRAYS

PRODUCT NO.	DESCRIPTION	PRODUCT NO.	DESCRIPTION
562-1100	Modular Tibial Tray Neutral Size 1	562-2100	Modular Tibial Tray Sloped Size 1
562-1200	Modular Tibial Tray Neutral Size 2	562-2200	Modular Tibial Tray Sloped Size 2
562-1300	Modular Tibial Tray Neutral Size 3	562-2300	Modular Tibial Tray Sloped Size 3
562-1400	Modular Tibial Tray Neutral Size 4	562-2400	Modular Tibial Tray Sloped Size 4
562-1500	Modular Tibial Tray Neutral Size 5	562-2500	Modular Tibial Tray Sloped Size 5
562-1600	Modular Tibial Tray Neutral Size 6	562-2600	Modular Tibial Tray Sloped Size 6

### TIBIAL INSERTS

PRODUCT NO.	DESCRIPTION	PRODUCT NO.	DESCRIPTION	PRODUCT NO.	DESCRIPTION
563-1108	CK Tibial Insert Size 1 8mm	563-1308	CK Tibial Insert Size 3 8mm	563-1508	CK Tibial Insert Size 5 8mm
563-1110	CK Tibial Insert Size 1 10mm	563-1310	CK Tibial Insert Size 3 10mm	563-1510	CK Tibial Insert Size 5 10mm
563-1112	CK Tibial Insert Size 1 12mm	563-1312	CK Tibial Insert Size 3 12mm	563-1512	CK Tibial Insert Size 5 12mm
563-1114	CK Tibial Insert Size 1 14mm	563-1314	CK Tibial Insert Size 3 14mm	563-1514	CK Tibial Insert Size 5 14mm
563-1116	CK Tibial Insert Size 1 16mm	563-1316	CK Tibial Insert Size 3 16mm	563-1516	CK Tibial Insert Size 5 16mm
563-1118	CK Tibial Insert Size 1 18mm	563-1318	CK Tibial Insert Size 3 18mm	563-1518	CK Tibial Insert Size 5 18mm
563-1120	CK Tibial Insert Size 1 20mm	563-1320	CK Tibial Insert Size 3 20mm	563-1520	CK Tibial Insert Size 5 20mm
563-1122	CK Tibial Insert Size 1 22mm	563-1322	CK Tibial Insert Size 3 22mm	563-1522	CK Tibial Insert Size 5 22mm
563-1124	CK Tibial Insert Size 1 24mm	563-1324	CK Tibial Insert Size 3 24mm	563-1524	CK Tibial Insert Size 5 24mm
563-1126	CK Tibial Insert Size 1 26mm	563-1326	CK Tibial Insert Size 3 26mm	563-1526	CK Tibial Insert Size 5 26mm
563-1128	CK Tibial Insert Size 1 28mm	563-1328	CK Tibial Insert Size 3 28mm	563-1528	CK Tibial Insert Size 5 28mm
563-1130	CK Tibial Insert Size 1 30mm	563-1330	CK Tibial Insert Size 3 30mm	563-1530	CK Tibial Insert Size 5 30mm
563-1208	CK Tibial Insert Size 2 8mm	563-1408	CK Tibial Insert Size 4 8mm	563-1608	CK Tibial Insert Size 6 8mm
563-1210	CK Tibial Insert Size 2 10mm	563-1410	CK Tibial Insert Size 4 10mm	563-1610	CK Tibial Insert Size 6 10mm
563-1212	CK Tibial Insert Size 2 12mm	563-1412	CK Tibial Insert Size 4 12mm	563-1612	CK Tibial Insert Size 6 12mm
563-1214	CK Tibial Insert Size 2 14mm	563-1414	CK Tibial Insert Size 4 14mm	563-1614	CK Tibial Insert Size 6 14mm
563-1216	CK Tibial Insert Size 2 16mm	563-1416	CK Tibial Insert Size 4 16mm	563-1616	CK Tibial Insert Size 6 16mm
563-1218	CK Tibial Insert Size 2 18mm	563-1418	CK Tibial Insert Size 4 18mm	563-1618	CK Tibial Insert Size 6 18mm
563-1220	CK Tibial Insert Size 2 20mm	563-1420	CK Tibial Insert Size 4 20mm	563-1620	CK Tibial Insert Size 6 20mm
563-1222	CK Tibial Insert Size 2 22mm	563-1422	CK Tibial Insert Size 4 22mm	563-1622	CK Tibial Insert Size 6 22mm
563-1224	CK Tibial Insert Size 2 24mm	563-1424	CK Tibial Insert Size 4 24mm	563-1624	CK Tibial Insert Size 6 24mm
563-1226	CK Tibial Insert Size 2 26mm	563-1426	CK Tibial Insert Size 4 26mm	563-1626	CK Tibial Insert Size 6 26mm
563-1228	CK Tibial Insert Size 2 28mm	563-1428	CK Tibial Insert Size 4 28mm	563-1628	CK Tibial Insert Size 6 28mm
563-1230	CK Tibial Insert Size 2 30mm	563-1430	CK Tibial Insert Size 4 30mm	563-1630	CK Tibial Insert Size 6 30mm

## CEMENTED STEMS

PRODUCT NO.	DESCRIPTION	PRODUCT NO.	DESCRIPTION
565-1203	Cemented Stem Extension 12mm x 30mm	565-1403	Cemented Stem Extension 14mm x 30mm
565-1205	Cemented Stem Extension 12mm x 50mm	565-1405	Cemented Stem Extension 14mm x 50mm
565-1208	Cemented Stem Extension 12mm x 80mm	565-1408	Cemented Stem Extension 14mm x 80mm
565-1212	Cemented Stem Extension 12mm x 120mm	565-1412	Cemented Stem Extension 14mm 120mm

## FLUTED STEMS

PRODUCT NO.	DESCRIPTION	PRODUCT NO.	DESCRIPTION
566-1008	Fluted Stem Extension 10mm x 80mm	566-1615	Fluted Stem Extension 16mm x 150mm
566-1011	Fluted Stem Extension 10mm x 110mm	566-1708	Fluted Stem Extension 17mm x 80mm
566-1015	Fluted Stem Extension 10mm x 150mm	566-1711	Fluted Stem Extension 17mm x 110mm
566-1208	Fluted Stem Extension 12mm x 80mm	566-1715	Fluted Stem Extension 17mm x 150mm
566-1211	Fluted Stem Extension 12mm x 110mm	566-1808	Fluted Stem Extension 18mm x 80mm
566-1215	Fluted Stem Extension 12mm x 150mm	566-1811	Fluted Stem Extension 18mm x 110mm
566-1308	Fluted Stem Extension 13mm x 80mm	566-1815	Fluted Stem Extension 18mm x 150mm
566-1311	Fluted Stem Extension 13mm x 110mm	566-2008	Fluted Stem Extension 20mm x 80mm
566-1315	Fluted Stem Extension 13mm x 150mm	566-2011	Fluted Stem Extension 20mm x 110mm
566-1408	Fluted Stem Extension 14mm x 80mm	566-2015	Fluted Stem Extension 20mm x 150mm
566-1411	Fluted Stem Extension 14mm x 110mm	566-2208	Fluted Stem Extension 22mm x 80mm
566-1415	Fluted Stem Extension 14mm x 150mm	566-2211	Fluted Stem Extension 22mm x 110mm
566-1508	Fluted Stem Extension 15mm x 80mm	566-2215	Fluted Stem Extension 22mm x 150mm
566-1511	Fluted Stem Extension 15mm x 110mm	566-2408	Fluted Stem Extension 24mm x 80mm
566-1515	Fluted Stem Extension 15mm x 150mm	566-2411	Fluted Stem Extension 24mm x 110mm
566-1608	Fluted Stem Extension 16mm x 80mm	566-2415	Fluted Stem Extension 24mm x 150mm
566-1611	Fluted Stem Extension 16mm x 110mm		

## TIBIAL AUGMENTS

PRODUCT NO.	DESCRIPTION	PRODUCT NO.	DESCRIPTION
567-1101	Tibial Augment 10° Wedge Rt Lat/Lt Med Size 1	567-2204	Tibial Augment 20° Wedge Lt Lat/Rt Med Size 4
567-1102	Tibial Augment 10° Wedge Rt Lat/Lt Med Size 2	567-2205	Tibial Augment 20° Wedge Lt Lat/Rt Med Size 5
567-1103	Tibial Augment 10° Wedge Rt Lat/Lt Med Size 3	567-2206	Tibial Augment 20° Wedge Lt Lat/Rt Med Size 6
567-1104	Tibial Augment 10° Wedge Rt Lat/Lt Med Size 4	567-5101	Tibial Augment 5mm Block Rt Lat/Lt Med Size 1
567-1105	Tibial Augment 10° Wedge Rt Lat/Lt Med Size 5	567-5102	Tibial Augment 5mm Block Rt Lat/Lt Med Size 2
567-1106	Tibial Augment 10° Wedge Rt Lat/Lt Med Size 6	567-5103	Tibial Augment 5mm Block Rt Lat/Lt Med Size 3
567-1201	Tibial Augment 10° Wedge Lt Lat/Rt Med Size 1	567-5104	Tibial Augment 5mm Block Rt Lat/Lt Med Size 4
567-1202	Tibial Augment 10° Wedge Lt Lat/Rt Med Size 2	567-5105	Tibial Augment 5mm Block Rt Lat/Rt Med Size 5
567-1203	Tibial Augment 10° Wedge Lt Lat/Rt Med Size 3	567-5106	Tibial Augment 5mm Block Rt Lat/Lt Med Size 6
567-1204	Tibial Augment 10° Wedge Lt Lat/Rt Med Size 4	567-5201	Tibial Augment 5mm Block Lt Lat/Rt Med Size 1
567-1205	Tibial Augment 10° Wedge Lt Lat/Rt Med Size 5	567-5202	Tibial Augment 5mm Block Lt Lat/Rt Med Size 2
567-1206	Tibial Augment 10° Wedge Lt Lat/Rt Med Size 6	567-5203	Tibial Augment 5mm Block Lt Lat/Rt Med Size 3
567-2101	Tibial Augment 20° Wedge Rt Lat/Lt Med Size 1	567-5204	Tibial Augment 5mm Block Lt Lat/Rt Med Size 4
567-2102	Tibial Augment 20° Wedge Rt Lat/Lt Med Size 2	567-5205	Tibial Augment 5mm Block Lt Lat/Rt Med Size 5
567-2103	Tibial Augment 20° Wedge Rt Lat/Lt Med Size 3	567-5206	Tibial Augment 5mm Block Lt Lat/Rt Med Size 6
567-2104	Tibial Augment 20° Wedge Rt Lat/Lt Med Size 4	567-6101	Tibial Augment 10mm Block Rt Lat/Lt Med Size 1
567-2105	Tibial Augment 20° Wedge Rt Lat/Lt Med Size 5	567-6102	Tibial Augment 10mm Block Rt Lat/Lt Med Size 2
567-2106	Tibial Augment 20° Wedge Rt Lat/Lt Med Size 6	567-6103	Tibial Augment 10mm Block Rt Lat/Lt Med Size 3
567-2201	Tibial Augment 20° Wedge Lt Lat/Rt Med Size 1	567-6104	Tibial Augment 10mm Block Rt Lat/Lt Med Size 4
567-2202	Tibial Augment 20° Wedge Lt Lat/Rt Med Size 2	567-6105	Tibial Augment 10mm Block Rt Lat/Lt Med Size 5
567-2203	Tibial Augment 20° Wedge Lt Lat/Rt Med Size 3	567-6106	Tibial Augment 10mm Block Rt Lat/Lt Med Size 6