

REDUCED-PROFILE
INTEGRATED® SPINE SYSTEM
NEW PRODUCT RELEASE

The Integrated® Spine System is designed for component strength and ease-of-use. This provides a design without compromising patient needs or surgeon satisfaction. The system includes the following features and benefits:



- Full 30° polyaxial screw design is versatile to rod placement and secure fixation
- Pre-bent rods in 5mm increments to reduce the need for intraoperative contouring
- Reduced-profile design and increased strength provide a sleeker fit and prevent splaying/cross-threading
- Minimal components required to complete instrumentation and implantation
- Persuasion up to 25mm
- Narrow head design for close screw placement
- Adaptable to other systems

ORDER INFORMATION

TRANSVERSE CONNECTOR

PRODUCT NO.	DESCRIPTION	PRODUCT NO.	DESCRIPTION
158-4050	Transverse Connector Short 40-50mm	158-6070	Transverse Connector Long 60-70mm
158-5060	Transverse Connector Medium 50-60mm	158-7080	Transverse Connector Extra Long 70-80mm

PEDICLE SCREW

PRODUCT NO.	DESCRIPTION	PRODUCT NO.	DESCRIPTION
158-5525	Pedicle Screw 5.5 x 25mm	158-6550	Pedicle Screw 6.5 x 50mm
158-5530	Pedicle Screw 5.5 x 30mm	158-6555	Pedicle Screw 6.5 x 55mm
158-5535	Pedicle Screw 5.5 x 35mm	158-6560	Pedicle Screw 6.5 x 60mm
158-5540	Pedicle Screw 5.5 x 40mm	158-7525	Pedicle Screw 7.5 x 25mm
158-5545	Pedicle Screw 5.5 x 45mm	158-7530	Pedicle Screw 7.5 x 30mm
158-5550	Pedicle Screw 5.5 x 50mm	158-7535	Pedicle Screw 7.5 x 35mm
158-6525	Pedicle Screw 6.5 x 25mm	158-7540	Pedicle Screw 7.5 x 40mm
158-6530	Pedicle Screw 6.5 x 30mm	158-7545	Pedicle Screw 7.5 x 45mm
158-6535	Pedicle Screw 6.5 x 35mm	158-7550	Pedicle Screw 7.5 x 50mm
158-6540	Pedicle Screw 6.5 x 40mm	158-7555	Pedicle Screw 7.5 x 55mm
158-6545	Pedicle Screw 6.5 x 45mm	158-7560	Pedicle Screw 7.5 x 60mm

CURVED SPINAL ROD

PRODUCT NO.	DESCRIPTION	PRODUCT NO.	DESCRIPTION
158-0030	Spinal Rod 5.5 x 30mm*	158-1030	Spinal Rod 6.0 x 30mm
158-0035	Spinal Rod 5.5 x 35mm*	158-1035	Spinal Rod 6.0 x 35mm
158-0040	Spinal Rod 5.5 x 40mm*	158-1040	Spinal Rod 6.0 x 40mm
158-0045	Spinal Rod 5.5 x 45mm*	158-1045	Spinal Rod 6.0 x 45mm
158-0050	Spinal Rod 5.5 x 50mm*	158-1050	Spinal Rod 6.0 x 50mm
158-0055	Spinal Rod 5.5 x 55mm*	158-1055	Spinal Rod 6.0 x 55mm
158-0060	Spinal Rod 5.5 x 60mm*	158-1060	Spinal Rod 6.0 x 60mm
158-0065	Spinal Rod 5.5 x 65mm*	158-1065	Spinal Rod 6.0 x 65mm
158-0070	Spinal Rod 5.5 x 70mm*	158-1070	Spinal Rod 6.0 x 70mm
158-0075	Spinal Rod 5.5 x 75mm*	158-1075	Spinal Rod 6.0 x 75mm
158-0080	Spinal Rod 5.5 x 80mm*	158-1080	Spinal Rod 6.0 x 80mm
158-0090	Spinal Rod 5.5 x 90mm*	158-1090	Spinal Rod 6.0 x 90mm
158-0105	Spinal Rod 5.5 x 105mm*	158-1105	Spinal Rod 6.0 x 105mm

* available upon request

STRAIGHT SPINAL ROD

PRODUCT NO.	DESCRIPTION	PRODUCT NO.	DESCRIPTION
158-6200	Spinal Rod 6.0 x 200mm	158-6405	Spinal Rod 6.0 x 405mm

CAP AND SET SCREW

PRODUCT NO.	DESCRIPTION
158-5000	Cap and Set Screw



12187 South Business Park Drive
Draper, Utah 84020

Telephone: 801.553.9991
Fax: 801.553.9993

For more information visit
www.orthodevelopment.com