



APPLICATION

Human musculoskeletal allograft may be used to repair bone defects, soft tissue injuries, and replacement of musculoskeletal tissue.

FORM

The Cervical Allograft Milled Cortical Wedge tissue has been retrieved under aseptic technique. The tissue is then processed aseptically and lyophilized (freeze-dried) for storage.

DESCRIPTION

The Cervical Allograft Milled Cortical Wedge is machined to a range of varying heights. See product description chart for height dimensions. All tissue has been collected, processed, stored and distributed in accordance with the American Association of Tissue Banks Standards (AATBS) and FDA Regulations. Please see the package insert for intended uses/indication, device description, contraindications, precautions, warnings and potential risks.

PACKAGING

The Cervical Allograft Milled Cortical Wedge is finally sterilized by exposure to gamma irradiation and supplied individual in a sterile peel package. Product labels are provided in quintuplet for convenience in record keeping.

To order call Ortho Development
Customer Service at
(800) 429-8339

Product No.	Description [†]
01-197-005	Milled Cortical Wedge 5mm
01-197-006	Milled Cortical Wedge 6mm
01-197-007	Milled Cortical Wedge 7mm
01-197-008	Milled Cortical Wedge 8mm
01-197-009	Milled Cortical Wedge 9mm
01-197-010	Milled Cortical Wedge 10mm
01-197-011*	Milled Cortical Wedge 11mm
01-197-012*	Milled Cortical Wedge 12mm
01-197-013*	Milled Cortical Wedge 13mm

* Sizes 11, 12 and 13mm are special orders.

† Actual height included in description.

Product Of:

