

Escalade® Acetabular System



## Escalade® Acetabular Shell

- True hemisphere design
- Titanium plasma spray coating
- Designed to maximize head size while maintaining bearing thickness
- Combination taper and snap locking mechanism
- Anti-rotational tabs are flush with the shell on all liners
- +0 liner face truly flush with the shell

## Legend® Polyethylene

- GUR 1020 resin utilized for superior mechanical properties
- 75kGy gamma irradiated crosslinking process
- Proprietary thermal treatment to eliminate free radicals
- EtO sterilized to maintain crosslinked integrity

## Femoral Head Compatibility

- Ortho Development CoCr heads: 28mm, 32mm, 36mm and 40mm
- BioloX® Delta Ceramic heads: 28mm, 32mm, 36mm and 40mm

## Legend® Neutral and Hooded Liners

		SHELL DIAMETER													
		44	46	48	50	52	54	56	58	60	62	64	66	68	70
HEAD DIAMETER	28	28/44	28/46	28/48	28/50	28/52	28/54	28/56	28/58	28/60	28/62	28/64-70			
	32		32/46	32/48	32/50	32/52	32/54	32/56	32/58	32/60	32/62	32/64-70			
	36				36/50	36/52	36/54	36/56	36/58	36/60	36/62	36/64-70			
	40						40/54	40/56	40/58	40/60	40/62	40/64-70			
HOODED LINERS AVAILABLE IN THESE SIZES								SIZES IN BOXES WITH BOLD OUTLINES ARE SPECIAL ORDER ONLY							

\*ALL NEUTRAL AND HOODED LINERS HAVE A 1 MM LATERALIZATION.

## Legend® +3 Neutral (Lateralized) and Face-Changing Liners

		SHELL DIAMETER													
		44	46	48	50	52	54	56	58	60	62	64	66	68	70
HEAD DIAMETER	28	28/44 +3	28/46 +3	28/48 +3	28/50 +3	28/52 +3	28/54 +3	28/56 +3	28/58 +3	28/60 +3	28/62 +3	28/64-70 +3			
	32		32/46 +3	32/48 +3	32/50 +3	32/52 +3	32/54 +3	32/56 +3	32/58 +3	32/60 +3	32/62 +3	32/64-70 +3			
	36				36/50 +3	36/52 +3	36/54 +3	36/56 +3	36/58 +3	36/60 +3	36/62 +3	36/64-70 +3			
	40						40/54 +3	40/56 +3	40/58 +3	40/60 +3	40/62 +3	40/64-70 +3			
								SIZES IN BOXES WITH BOLD OUTLINES ARE SPECIAL ORDER ONLY							

\*LEGEND AND ESCALADE LINERS ARE COMPATIBLE WITH ALL ORTHO DEVELOPMENT SHELLS.