

BKS TriMax[®]



 **ORTHO**
DEVELOPMENT[®]

BKS TriMax® Reference Guide




Tibial Inserts are available in CR, UC, and PS in 6 sizes with thicknesses of 7mm, 8mm, 9mm, 10mm, 11mm, 12mm, 13mm, 14mm and 16mm (18mm and 20mm UC and PS only).

BKS TriMax® Femoral Components are compatible with E-Vitalize® Tibial Inserts only.

Tibial Inserts are interchangeable with the Femoral Components and Tibial Trays, see chart below.

| | | FEMORAL COMPONENT SIZE | | | | | |
|------------------|---|------------------------|---|------|------|---|---|
| | | 2 | 3 | 4N/4 | 5N/5 | 6 | 7 |
| TIBIAL TRAY SIZE | 2 | 2 | 2 | | | | |
| | 3 | 3 | 3 | 3 | | | |
| | 4 | | 4 | 4 | 4 | | |
| | 5 | | | 5 | 5 | 5 | |
| | 6 | | | | 6 | 6 | 6 |
| | 7 | | | | | 7 | 7 |

Shaded areas indicate Tibial Insert sizes.



FEMORAL COMPONENT

| SIZE | A/P (MM) | M/L (MM) |
|------|----------|----------|
| 2 | 54.0 | 59.5 |
| 3 | 57.5 | 62.5 |
| 4N | 61.5 | 62.5 |
| 4 | 61.5 | 66.5 |
| 5N | 65.5 | 66.5 |
| 5 | 65.5 | 70.5 |
| 6 | 69.5 | 74.5 |
| 7 | 74.0 | 79.5 |

TIBIAL TRAY

| SIZE | A/P (MM) | M/L (MM) |
|------|----------|----------|
| 2 | 39.0 | 61.0 |
| 3 | 41.5 | 65.0 |
| 4 | 44.5 | 69.5 |
| 5 | 48.0 | 75.0 |
| 6 | 51.5 | 80.5 |
| 7 | 55.0 | 86.0 |

PATELLA

| SIZE (DIAMETER) | HEIGHT (THICKNESS) |
|--------------------|-----------------------|
| 29mm | 7.5mm |
| 32mm | 8.0mm |
| 35mm | 9.0mm |
| 38mm | 10.0mm |
| 41mm | 11.0mm |



For more information about the BKS TriMax[®] please contact Ortho Development.

Please see the Balanced Knee[®] System Instructions for Use for intended use/indications, device description, contraindications, precautions, warnings, and potential risks associated with the BKS TriMax[®].



Ortho Development Corporation
odev.com

12187 So. Business Park Drive
Draper, Utah 84020
801-553-9991/fax 801-553-9991

Copyright © 2016 Ortho Development[®] Corporation. All rights reserved.
Unless otherwise indicated, all materials on these pages are copyrighted by Ortho Development Corporation. All rights reserved. No part of these pages, either text or image may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, electronic, mechanical or otherwise, for reasons other than personal use, is strictly prohibited without prior written permission.