




# Envision®

The Envision® Anterior Cervical Plate offers a unique blend of features engineered for ease-of-use and superior clinical results. The rotationally-dynamic design enables load-sharing between implant and graft for a more reliable fusion. It has a seamless, one-step self locking mechanism while maintaining a low profile (2.0mm) that reduces soft tissue irritation.

- Patented one-step locking mechanism employs both a rim and a taper lock
- Large graft windows provide superior visualization of the interbody space
- Screw allows for variable angulation up to 15° in any direction (30° total)

## PLATES

ONE LEVEL	TWO LEVEL	THREE LEVEL	FOUR LEVEL
14mm	26mm	46mm	60mm
16mm	28mm	49mm	65mm
18mm	30mm	52mm	70mm
20mm	32mm	55mm	75mm
22mm	35mm	58mm	80mm
24mm	38mm	61mm	85mm
	41mm	64mm	90mm
	44mm	67mm	95mm
	47mm		



## SCREWS

	4.0 MM (DIA.)	4.35 MM (DIA.)
12mm	Green	Dark Green
14mm	Blue	Dark Blue
16mm	Magenta	Dark Magenta
18mm*	Bronze	Dark Bronze

\*available by special order only

