

# Balanced Knee® System



 **ORTHO**  
DEVELOPMENT®

## Balanced Knee® System Reference Guide

Tibial Inserts are available in PS, CR and UC

7 sizes with thicknesses of 7mm, 8mm, 9mm, 10mm, 11mm, 12mm, 13mm, 14mm, 16mm, 18mm and 20mm





Tibial Inserts are interchangeable with the Femoral Components and Tibial Trays, see chart below

- Matched one-to-one with the Tibial Tray
- Up/Down one size with the Femoral Component

### FEMORAL COMPONENT SIZE

		1	2	3	4	4	5	5	6	7
					Narrow		Narrow			
TIBIAL TRAY SIZE	1	1	1							
	2	2	2	2						
	3		3	3	3	3				
	4			4	4	4	4	4		
	5				5	5	5	5	5	
	6						6	6	6	6
	7								7	7

Shaded areas indicate Tibial Insert sizes.



### FEMORAL COMPONENT

SIZE	A/P (MM)	M/L (MM)
1	50.0	56.5
2	54.0	59.5
3	57.5	62.5
4N	61.5	62.5
4	61.5	66.5
5N	65.5	66.5
5	65.5	70.5
6	69.5	74.5
7	74.0	79.5

### TIBIAL TRAY

SIZE	A/P (MM)	M/L (MM)
1	36.5	57.5
2	39.0	61.0
3	41.5	65.0
4	44.5	69.5
5	48.0	75.0
6	51.5	80.5
7	55.0	86.0

### PATELLA

SIZE (DIAMETER)	HEIGHT (THICKNESS)
29mm	7.5mm
32mm	8.0mm
35mm	9.0mm
38mm	10.0mm

For more information about the Balanced Knee® System please contact Ortho Development.

Please see the Balanced Knee® System Surgical Technique and Instructions for Use for intended use/indications, device description, contraindications, precautions, warnings, and potential risks associated with the Balanced Knee® System.



Ortho Development Corporation  
odev.com

12187 So. Business Park Drive  
Draper, Utah 84020  
801-553-9991/fax 801-553-9991

Copyright © 2014 Ortho Development® Corporation. All rights reserved.

Unless otherwise indicated, all materials on these pages are copyrighted by Ortho Development Corporation. All rights reserved. No part of these pages, either text or image may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, electronic, mechanical or otherwise, for reasons other than personal use, is strictly prohibited without prior written permission.